



MedSource
Diagnostics

Electrodes



The MedSource Foam Electrodes use a conductive adhesive hydro-gel “wet gel” which allows for maximum adhesion and electrical contact.

Our electrodes are made of fluid and tear resistant foam. Special tabs help for easy removal and application.

The foam electrodes do not contain natural rubber latex.



www.gomedsource.com

Rev. 1

Item #	Description	UOM
MS-65003	Electrodes, 3/PK, 600 EA/CS	CS
MS-65005	Electrodes, 30/PK, 600 EA/CS	CS