



**V-VAC Starter Kit for Patient Use (shown above) includes:**  
 Handle, 2 Replacement Cartridges, Short Suction Catheter w/Adapter Tip, Double Male Connector, and Directions for Use.

Ordering information:

- 98 50 00 V-VAC Starter Kit
- 98 51 00 V-VAC Training Kit
- 98 50 01 Replacement Cartridge
- 98 52 00 Carrying Bag
- 98 50 30 Replacement Handle

**V-VAC Training Kit for Practical Use includes:**

Handle, 1 Training Cartridge with plug for emptying and cleaning, Catheter w/Adapter Tip and Double Male Connector.\*

\*Not illustrated.



**A** The V-VAC unit fits easily into many emergency kits. A separate carry bag is also available.

**B** The soft tip of the intake valve has deep cutouts to prevent tube occlusion and soft tissue injury. The tri-cut valve resists back flow of suctioned matter.

For additional information or to order, contact your local Laerdal Company or Distributor



**Norway:** LAERDAL MEDICAL AS  
 PO Box 377, N-4002 Stavanger  
 Tel: +47 51 51 17 00, Fax +47 51 52 35 57  
 E-mail: laerdal.norway@laerdal.no

**Australia:** LAERDAL Pty Ltd. (ABN 47 003 817 490)  
 23 Edward Street, Oakleigh, Victoria 3166  
 Tel: +61 3 9569 4055, Toll free 1800 331 565  
 Fax +61 3 9569 4028, Toll free 1800 635 835  
 E-mail: customerservice@laerdal.com.au

**Benelux:** "The Laerdal BeNeLux House"  
 Dijkstraat 54, NL-5554 PS VALKENSWAARD  
 Tel: +31 40 208 58 00, Fax +31 40 208 58 02  
 E-mail: info@laerdal.nl

**Benelux:** Belgium and Luxembourg  
 LAERDAL BeNeLux nv  
 Mechelsesteenweg 277, B-1800 VILVOORDE  
 Tel: +32 2 253 36 96, Fax + 32 2 253 36 90  
 E-mail: info@laerdal.be

**Canada and Latin America:**  
 LAERDAL MEDICAL CANADA LTD.  
 151 Nashdene Rd., Unit #45  
 Toronto, ON, Canada, M1V 4C3  
 Tel: +1 (416) 298-9600  
 Toll free 888/LAERDAL (523-7325)  
 ou en français (800) 567-9987  
 Fax +1 (416) 298-8016  
 E-mail: savelves@laerdal.ca

**Denmark:** LAERDAL DANMARK  
 Esplanaden 18, 1  
 DK-1263 Copenhagen K  
 Tel: +45 80 333 112, Fax: +45 80 333 555  
 E-mail: laerdal.denmark@laerdal.no

**Eastern Europe:** LAERDAL C.E. EUROPE  
 A-1060 Vienna, Mariahilferstr. 1/d/1/3  
 Tel-Fax +43 15877140, Tel: +43 15815927  
 E-mail: mstam@aon.at

**Far East:** LAERDAL Singapore Pte Ltd  
 No 1 Marine Parade Central  
 #13-05, Parkway Centre, Singapore 449408  
 Tel: +65 63464259, Fax: +65 63467523  
 E-mail: admin@laerdal.com.sg

**Finland:** LAERDAL OY  
 Vattunienranta 2, 00210 HELSINKI  
 Puhelin +358 (0)9-612 99 80,  
 Fax +358 (0)9-692 77 99  
 E-mail: laerdal.finland@laerdal.no

**France:** LAERDAL MEDICAL FRANCE  
 1 rue des Vergers - Bâtiment n° 5 - 69578 Limonest Cedex  
 Tel: +33 (0)4 72 52 02 52, Fax: +33 (0)4 78 35 38 45  
 E-mail: laerdal.france@laerdal.no

**Germany:** LAERDAL MEDICAL DEUTSCHLAND GmbH  
 Am Lohrlefeld 56, 81249 München  
 Tel: +49 (0)89 / 864 954-0, Fax +49 (0)89 / 864 34 84  
 E-mail: info@laerdal.de

**Italy:** LAERDAL ITALIA s.r.l.  
 Via Piero Gobetti, 5/22, int.Z, 40129 Bologna  
 Tel: +39 051-355587, Fax: +39 051-355598  
 E-mail: laerdal.italy@laerdal.no

**Japan:** LAERDAL MEDICAL JAPAN K.K.  
 Shinjuku IS Bldg., 8F  
 2-13-12, Shinjuku, Shinjuku-ku, Tokyo 160-0022  
 Tel: 81-3-5919-2121, Fax: 81-3-5919-2122  
 E-mail: laerdal.japan@laerdal.co.jp

**Malaysia:** LAERDAL HOSPILINE Sdn Bhd  
 12 Jalan Titar 33/25, Section 33, Shah Alam Technology Park,  
 40400 Shah Alam, Selangor Darul Ehsan,  
 Tel: +60 (3) 5122-7002, Fax +60 (3) 5122-5650  
 E-mail: customerservice@laerdal.com.my

**New Zealand:** LAERDAL NEW ZEALAND Ltd  
 (GST Number 79-403-636)  
 Po Box 36339, Merivale, Christchurch  
 Tel: +61 3 9564 4066, Toll free 0800 523 732  
 Fax +61 3 9563 3368, Toll free 0800 528 852  
 E-mail: customerservice@laerdal.co.nz

**Portugal:** LAERDAL ESPANA, Sucursal em Portugal  
 Rua da Bela Vista à Graca, no.31 A, Escritorio 9,  
 1170-054 Lisboa  
 Tel: (+351)-21-816-6567, Fax (+351)-21-816-6568,  
 Mob. (+351)-918765367

**Spain:** LAERDAL ESPANA, S.L.  
 Manuel Tovar 19, 28034 Madrid  
 Tel: (+34) 902-291110, Fax: (+34) (91)-7291125  
 E-mail: laerdal.spain@laerdal.no

**Sweden:** LAERDAL MEDICAL AB  
 Box 2064, Paviljongvägen 5, 132 02 Saltsjö-Boo  
 Tel: +46 (0)8-55614610, Fax +46 (0)8-55614619  
 E-mail: kontakt@laerdal.se

**United Kingdom:** LAERDAL MEDICAL LTD.  
 Laerdal House, Goodmead Road, Orpington,  
 Kent BR6 0HX  
 Tel: +44 (0)1689 876634, Fax +44 (0)1689 873800  
 E-mail: customerservice@laerdal.co.uk

**USA:**  
 LAERDAL MEDICAL CORPORATION  
 167 Myers Corners Road, PO Box 1840  
 Wappingers Falls, New York 12590-8840  
 Tel: (800) 431-1055, +1 (845) 297-7770  
 Fax: (800) 227-1143, +1 (845) 298-4545  
 E-mail: customerservice@laerdal.com

**LAERDAL TEXAS**  
 PO Box 38, 226 FM 116 South Gatesville,  
 Texas 76528-0038  
 Tel: (800) 433-5539, (254) 865-7221  
 Fax: (254) 865-8011  
 E-mail: customerservice@laerdal.com

Get a good grip on handling even the toughest suction scenario



## V-VAC Manual Suction Unit

Ideal for First Responder emergency kits and as back-up suction for hospitals, clinics and emergency responders

- All-in-one system adjusts for high and low vacuum suctioning
- Efficient: exceeds 380 mmHg suction with a peak air flow of over 70 litres per minute
- Unique, soft, no clog suction tip
- Catheter connection allows for versatile use
- Disposable, one-piece cartridge eliminates the need for cleaning contaminated parts
- Reusable handle is easy to clean for repeated use

